



## ULTRA CLEANING POWER FOR YOUR CSSD

### Aniosyme Synergy WD Ultra

#### DI-ENZYMATIC DETERGENT FOR ULTRA CLEANING POWER

Aniosyme Synergy WD Ultra is the first combined ultra concentrate detergent for automated Instrument Reprocessing. The unique combined formula allows an adjustable pH-level from neutral to alkaline. Its di-enzymatic, silicate free formula offers excellent material compatibility with its unique low concentration.

- Combined neutral to alkaline pH-level formula
- Excellent cleaning performance starting from 0.03%
- Suitable above pH10 with any water quality
- Space saving - no more big canisters



# Aniosyme Synergy WD Ultra

Di-enzymatic detergent for alkaline cleaning power



Aniosyme Synergy WD Ultra  
5 L Canister



Aniosyme Synergy WD Ultra  
20 L Canister

## A UNIQUE FORMULATION

- Combined neutral to alkaline pH-level formula
- Excellent cleaning performance starting from 0.03%
- Suitable above pH10 with any water quality

## CHEMO-PHYSICAL DATA

- |                                                             |      |            |
|-------------------------------------------------------------|------|------------|
| • Density                                                   | 20°C | 1,158 g/mL |
| • pH value* (concentrate)                                   | 20°C | 9,9        |
| • pH value* at 0.05% (solution, 0.5 ml/l in DI water)       | 20°C | 10,2       |
| • pH value* at 1% (solution, 10 ml/l in DI water)           | 20°C | 11,0       |
| • Refractive index (concentrate)                            | 20°C | 1,4022     |
| • Conductivity at 0.05%<br>(solution, 0.5 ml/l in DI water) | 20°C | 85 µS/cm   |
| • Conductivity at 1% (solution, 10 ml/l in DI water)        | 20°C | 1421 µS/cm |
| • Viscosity (concentrate)                                   | 20°C | < 20 mPa/s |

\* pH value measurements under operating conditions always by means of a pH meter with temperature compensation.

## ORDERING INFORMATION

| PRODUCT                            | PACKAGING UNIT | ORDER CODE |
|------------------------------------|----------------|------------|
| Aniosyme Synergy WD Ultra Canister | 2 x 5 L        | XXXXXXX    |
| Aniosyme Synergy WD Ultra Canister | 1 x 20 L       | XXXXXXX    |



**ECOLAB DEUTSCHLAND GmbH**  
Ecolab-Allee 1  
40789 Monheim am Rhein  
Germany

Distributed by:  
**ECOLAB DEUTSCHLAND GmbH**  
Ecolab-Allee 1  
40789 Monheim am Rhein  
Germany

## APPLICATION AREA

Concentrated di-enzymatic detergent for instrument reprocessing. Cleaning of surgical and medical equipment, medical devices in automatic washing disinfectant and washing tunnel.

## INSTRUCTIONS FOR USE

- Concentrated solution for professional use.
- Use in an automatic washing machine and washing tunnel at an application dose of 0,05%. If necessary, the dosing can be adjusted between 0,03% to 1%.
- Effective from 1 minute.
- Adapt the concentration, the contact time and the temperature of the dilution water according to the degree of soiling, the cleaning method used and the targeted pH-level of the cleaning phase.
- Maximum temperature: 60°C.

## INGREDIENTS

Non-ionic surfactants, sequestering agent, enzyme complex (protease, lipase), excipients.

## STORAGE


Dangerous - respect the precautions for use (drawn up according to the European rules in force regarding the classification and labelling of chemical products). Storage: between +5°C and +25°C.

## CERTIFICATES

On demand we are pleased to provide you with certificates.

## MEDICAL DEVICE USE SAFELY

Please read label and product information.



FOR PROFESSIONAL USE ONLY.

