

3-COMPARTMENT SINK



SIMPLE. SAFE. SUSTAINABLE.

SolidSense™

- ▲ Powerful grease-cutting in a safe and easy to handle solid block
- ▲ Broad range EPA-registered solid sanitizer approved from 150–400 ppm
- ▲ Shape and color-coded blocks to maximize ease of use
- ▲ Front-loading dispenser for easy refill
- ▲ 64% less tonnage lowers shipping and storage costs
- ▲ 67% less packaging waste
- ▲ Patented dispensing technology ensures consistent dilution

A flexible, all solid 3-compartment sink program delivers performance, ease of use, and savings retailers require.



SolidSense™

Leveraging proven solid products with the latest advancements in solid dispensing technology, SolidSense now provides **improved flexibility, ease of use, and cost controls** that are critical for food retailers today.


Patented Solid Dispensing Technology




Solid block format

- Allows magazine loading to prevent running out of product.
- Improves dilution consistency by maintaining a constant distance to the water source.

Pressure and temperature regulators further reduce the impact of pressure and temperature changes, improving dilution consistency.

Turbulent flow technology mitigates impact of water pressure changes, improves dilution consistency.



PRODUCT	ITEM #	PACK SIZE
 SolidSense™ All Purpose Super Concentrate	1114170	4 x 2 lb blocks
 SolidSense™ Low Foam Powersink Detergent	1110195	4 x 2.45 lb blocks
 KAY® SolidSense™ Sanitizer	1114169	4 x 2 lb blocks

3-COMPARTMENT SINK

SIMPLE. SAFE. SUSTAINABLE.

BENEFITS

67% less packaging waste.



1 lb packaging waste = 3 lbs packaging waste

An average supermarket will save +80 lbs of packaging waste annually

BENEFITS

64% less tonnage shipped.



2.8 lbs liquid product

1 lb solid product

6144 wash sinks per pallet
1771.2 lbs

9402 wash sinks per pallet
1045.5 lbs

4915 sanitizer sinks per pallet
1704 lbs

9874 sanitizer sinks per pallet
961.2 lbs

An average supermarket will ship +1,200 fewer lbs annually

BENEFITS

39% fewer pallets stored.



Every pallet of solid product delivers:
3,200 more sinks of pot and pan detergent
5,000 more sinks of sanitizer

Savings estimates based on average annual pot and pan detergent and sanitizer usage of a typical supermarket. Sinks per pallet estimates based on filling sinks with 20 gal. of solution.

Kay Chemical Company
Food Retail Services
8300 Capital Drive
Greensboro, NC 27409-9790
800.529.5458

45846-00 FRS 45846/8000/0114 ©2014 Kay Chemical Company. All rights reserved.

